


Now FDA-Authorized

### Individual Patient Results

At Month 13, a significant reduction in drusen volume and no visible loss of photoreceptor/retinal pigment epithelium cells were observed.

Age: 77 years      Baseline BCVA\*: 75 letters  
 Sex: Female      Month 13 BCVA: 79 letters  
                          Month 21 BCVA: 84 letters  
                          Month 24 BCVA: 82 letters



Individual patient results may vary

\*OCT imaging and BCVA measurement taken at screening visit

## Valeda Demonstrates Improvements in Clinical and Anatomical Outcomes Supporting a Disease-Modifying Benefit

- PBM demonstrates a benefit in BCVA versus Sham over the course of the trial. The primary BCVA endpoint at Month 21 had a p value = 0.0036
- More subjects lost BCVA in the Sham group compared to the PBM group at Months 13, 21, and 24
- A greater numerical increase in macular drusen volume was observed in the Sham group versus the PBM group
- Incident GA was observed in 24.0% of Sham versus 6.8% of PBM-treated eyes at Month 24†
- A favourable safety profile was observed with no signs of phototoxicity

### Indications for Use

The indicated use is for treatment of ocular damage and disease using photobiomodulation, including inhibition of inflammatory mediators, edema or drusen deposition, improvement of wound healing following ocular trauma or surgery, and increase in visual acuity and contrast sensitivity in patients with degenerative diseases such as dry age-related macular degeneration.



# The First Approved Treatment for Dry AMD

It's Time for Patients to See Their Future



+44-207-660-1516  
 myvaleda.com



©2025 LumiThera, Inc.  
 MKT-0136 enGB REV A

# LIGHTSITE III

Double-masked, randomized, sham-controlled, parallel group, multi-center trial to assess the safety and efficacy of photobiomodulation (PBM) treatment with Valeda in subjects with dry age-related macular degeneration (AMD)

## Baseline Characteristics



**Subjects:** 100 (98 subjects mITT analysis)  
**Eyes:** 148 (145 eyes mITT analysis)  
**Randomization:** 2:1 PBM to Sham  
**Race:** 99% Caucasian; 1% Black/African American  
**Gender:** 32 Males (32%); 68 Females (68%)  
**Mean Age:** 75 years  
**Mean Time from Diagnosis:** 4.9 years  
**AREDS Supplements:** 86 (86%) yes; 14 (14%) no  
**BCVA Baseline (BL)  $\geq 70$  Letters (20/40):** 103 eyes (70%)  
**BCVA Letter Score:** PBM: 70.7 letters (SD 5.2); Sham: 70.1 letters (SD 4.3)

## LIGHTSITE III Trial Design

**PBM:** 590, 660, and 850 nm wavelengths  
**Sham:** 10x reduction of 590 nm, 100x reduction of 660 nm, and no 850 nm wavelengths


Starting BCVA between 20/32 - 20/100						Month 13 Analysis <sup>1</sup>	Month 21 Final Tx Visit Analysis <sup>1</sup>	Month 24 Final Visit
Tx Series 1	Tx Series 2	Tx Series 3	Tx Series 4	Tx Series 5	Tx Series 6			
PBM	PBM	PBM	PBM	PBM	PBM			
Sham	Sham	Sham	Sham	Sham	Sham			
9 Tx Sessions/3-5 Weeks	9 Tx Sessions/3-5 Weeks	9 Tx Sessions/3-5 Weeks	9 Tx Sessions/3-5 Weeks	9 Tx Sessions/3-5 Weeks	9 Tx Sessions/3-5 Weeks			

<sup>1</sup>Co-primary endpoints: 13- and 21-Month comparison between PBM and Sham groups. This trial summary includes data from Month 13, Month 21, and Month 24 (3 months following last treatment).

## Valeda Improves and Maintains Vision

- PBM demonstrates a benefit in BCVA versus Sham over the course of the trial. The primary BCVA endpoint at Month 21 had a p value = 0.0036
- PBM improves BCVA with a mean 6.0 letter gain from BL at Month 13 (p < 0.0001) and maintains a mean 6.2 letter gain from BL at Month 21 (p < 0.0001)


### BCVA Letter Gain



80 subjects/113 eyes completed through Month 24


The least squares (LS) mean data using multiple imputation and standard error (SE) are presented. \* p < 0.05 between group comparison; ^ p < 0.0001 within group comparison (PBM)

### BCVA Letter Gain Distribution at Month 21\*



\*Data presented with multiple imputation at Month 21


## BCVA >5 Letter Loss Over 24 Months\*




\*Number of actual eyes used for percentage of >5 BCVA letter loss at Month 24

## A Greater Numerical Increase in Macular Drusen Volume Observed in Sham Group Versus PBM Group

\*Data presented with multiple imputation at Month 21



## Incident Geographic Atrophy (GA) was Higher in Sham Group Versus PBM Group at Months 13 and 24<sup>†</sup>



<sup>†</sup>Incident GA was not a pre-specified endpoint. Month 13 (p = 0.024, Fisher exact test, odds ratio 9.4) and Month 24 (p = 0.007, Fisher exact test, odds ratio 4.2)